



Patella Band™

Single Patient Use Only



IMPORTANT

Read carefully and completely before applying.

INTENDED USE

This patella tendonitis band is a unique device that provides continuous compression and support to the patellar tendon.

INDICATIONS

Patellar Tendonitis; Jumpers' Knee and Osgood-Schlatter's Disease

PRECAUTIONS

Consult your medical professional immediately if you experience sensation changes, unusual reactions, swelling or prolonged pain while using this product.

CONTRAINDICATIONS

The Patella Band should not be used by patients with poor peripheral circulation, diabetes or decreased skin sensitivity.

Sizing: Measure knee circumference below kneecap. If you measure between sizes, order the smaller size.

X-Small	11" – 12"
Small	12" – 13"
Regular	13" – 14"
Large	14" – 15"
X-Large	15" – 16"
XX-Large	16" – 17"

DIRECTIONS FOR USE

- Stand or sit with leg slightly bent and thigh relaxed. Slip Patella Band onto leg placing the tubular buttress around






the knee joint just below the kneecap. Making sure the contoured part of the band fits comfortably just below the popliteal region (back of knee).

- Feed each strap through its appropriate D-Ring closure.
- To tighten, pull both straps simultaneously and attach to the band for a snug but comfortable fit. Trim excess strapping where needed.

WASHING INSTRUCTIONS

Hand wash in lukewarm water with mild detergent. Air Dry. DO NOT PUT IN THE DRYER.



Manufacturer
DeRoyal
200 DeBusk Lane, Powell, TN 37849 U.S.A.
888.938.7828

Authorized Representative
DeRoyal Europe, Ltd.
Virginia Road, Kells, County Meath Ireland
353.46.9280100

©2003 DeRoyal. All Rights Reserved.
Part # 0-1414 Revision 2 8/03



Patella Band™

Single Patient Use Only



IMPORTANT

Read carefully and completely before applying.

INTENDED USE

This patella tendonitis band is a unique device that provides continuous compression and support to the patellar tendon.

INDICATIONS

Patellar Tendonitis; Jumpers' Knee and Osgood-Schlatter's Disease

PRECAUTIONS

Consult your medical professional immediately if you experience sensation changes, unusual reactions, swelling or prolonged pain while using this product.

CONTRAINDICATIONS


The Patella Band should not be used by patients with poor peripheral circulation, diabetes or decreased skin sensitivity.

Sizing: Measure knee circumference below kneecap. If you measure between sizes, order the smaller size.

X-Small	11" – 12"
Small	12" – 13"
Regular	13" – 14"
Large	14" – 15"
X-Large	15" – 16"
XX-Large	16" – 17"

DIRECTIONS FOR USE

- Stand or sit with leg slightly bent and thigh relaxed. Slip Patella Band onto leg placing the tubular buttress around






the knee joint just below the kneecap. Making sure the contoured part of the band fits comfortably just below the popliteal region (back of knee).

- Feed each strap through its appropriate D-Ring closure.
- To tighten, pull both straps simultaneously and attach to the band for a snug but comfortable fit. Trim excess strapping where needed.

WASHING INSTRUCTIONS

Hand wash in lukewarm water with mild detergent. Air Dry. DO NOT PUT IN THE DRYER.



Manufacturer
DeRoyal
200 DeBusk Lane, Powell, TN 37849 U.S.A.
888.938.7828

Authorized Representative
DeRoyal Europe, Ltd.
Virginia Road, Kells, County Meath Ireland
353.46.9280100

©2003 DeRoyal. All Rights Reserved.
Part # 0-1414 Revision 2 8/03