EN: INSUFFLATION TUBING

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FR: TUYAU D'INSUFFLATION

DE: INSUFFLATIONSSCHLAUCH

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NL: INSUFFLATIESLANG

SE: UPPBLÅSNINGS SLANG

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MANUFACTURED BY:
DeRoyal Industries, Inc.
200 DeBusk Lane
Powell, TN 37849 USA
888.938.7828 or (001) 865.938.7828
www.deroyal.com

MedEnvoy Global B.V.
Prinses Margrietplantsoen 33 – Suite 123
2595 AM the Hague
The Netherlands

PART#0-1876G | REVISED 10/2023

©2023 DEROYAL INDUSTRIES, INC.
All Rights Reserved. DeRoyal and the DeRoyal logo are trademarks or registered trademarks of DeRoyal Industries, Inc.

EN: ENGLISH

INSUFFLATION TUBING

	MEDICAL DEVICE
	DO NOT REUSE
	STERILIZED WITH ETHYLENE OXIDE
	NOT MADE WITH NATURAL RUBBER LATEX
	FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.

IMPORTANT INFORMATION

Please read all instructions, warnings, and cautions before use. Correct application is essential for proper functioning of product.

INTENDED USE

DeRoyal® Insufflation Tubing is intended for use by a licensed practitioner to deliver carbon dioxide gas from the insufflator to the patient to assist in creating pneumoperitoneum during a laparoscopic procedure.

WARNINGS

- Insufflation tubing is not intended to be used for introduction of fluid into the body.
- This system is only recommended for laparoscopic surgery.
- Tubing sets will be connected to a trocar or veress needle and should not be used or inserted directly into the body, organs, and tissues.
- Inspect the product and packaging for visible cracks or broken components.
- DO NOT reuse or attempt to re-sterilize.
- Usage of this device other than directed by a physician or in these instructions may result in device malfunction and/or patient injury.
- DO NOT use if device causes patient to experience skin irritation.
- DO NOT use if filter or housing is damaged.
- DO NOT use above recommended insufflator settings.
- Purge device of air prior to delivering gas to patient.

CAUTIONS

This device is for single patient use only. It is not intended for disinfection and subsequent reuse, which may result in device failure and/or create the risk of contamination.

DIRECTIONS

- Open the insufflation tubing set in the sterile field using aseptic technique.
- The adapter or luer lock closest to the filter should be connected to the insufflator.
- This connection should be made with the luer lock or an adapter that is compatible with the insufflator.
- The opposite end of the tubing should be connected to either a trocar or veress needle (pneumo needle) for initial set-up. This

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connection should be made with the luer lock or an appropriate adapter.

- Immediately before delivering gas to the patient, purge air from tubing sets and instruments with CO2.
- If using a veress needle, move the insufflation tubing from the veress needle to the fitting on the trocar after pneumoperitoneum is established.

STERILIZATION INFORMATION

DeRoyal intends that non-sterile products for use in sterile environments be further processed by our customers, including further packaging and/or sterilization. These products are to be sterilized by ethylene oxide, utilizing the following parameters: Temperature: 120°F - 135°F (48.9°C - 57.2°C), Ethylene Oxide Gas Concentration: 625 - 800 mg/l, Exposure Time: 4 hrs minimum (Preconditioning: Temperature 100°F - 120°F (37.8°C - 48.9°C), Humidity 55-75%, Time 18hrs minimum), (Aeration: Temperature 100°F - 135°F (37.8°C - 42.2°C), Time 18hrs minimum). Other sterilization methods have not been validated by DeRoyal and may damage the product which could cause a device malfunction and/or injury to the patient.

The device has not been evaluated for re-processing or re-sterilization. Reprocessing and/or re-sterilization may damage the device, rendering it unusable and/or may lead to device failure which could result in patient illness, injury or death.

Product provided in sterile and non-sterile configurations. This is identified on the packaging with the following symbols:

	NON-STERILE
	STERILIZED WITH ETHYLENE OXIDE

STORAGE AND TRANSPORT CONDITIONS

Do not store product at extreme temperatures or in a moist/damp environment.

	KEEP DRY
	KEEP AWAY FROM SUNLIGHT

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

WARRANTY

DeRoyal® products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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