




## SILVER ALGINATE CATHETER DRESSING

# ALGIDEX™ Ag I.V. PATCH

### ALGIDEX Ag I.V. PATCH



### FEATURES & BENEFITS

- Reduces bacterial colonization associated with catheter related infections
- Absorbs moisture around the catheter insertion site
- Sterile patch of polyurethane foam coated with an ionic silver alginate and maltodextrin matrix
- Releases active ionic silver for broad antimicrobial effectiveness

### INDICATIONS


- Peripheral IV catheters
- Central venous lines
- Arterial catheters
- Dialysis catheters
- External fixator pins
- Peritoneal dialysis catheters
- Epidural catheters
- Gastrostomy feeding tubes

TO PLACE AN ORDER OR FOR MORE INFORMATION  
PLEASE CALL: **800-251-9864**  
865.938.7828 (OUTSIDE USA)  
800.327.0304 (FAX)

### APPLICATION & REMOVAL


#### APPLICATION:

1. May be applied immediately following initial catheter placement or during routine catheter dressing change.
2. Apply Algidex Ag I.V. Patch with dark side touching the catheter and skin.
3. Secure patch in place with retention dressing such as gauze, transparent film or tape.
4. Patch may be worn up to 7 days.



#### EXPECTED PERFORMANCE

When moisture or drainage is present at the catheter site, the Algidex matrix will separate from the foam disc as it absorbs the moisture/drainage.



#### REMOVAL:

1. Remove retention dressing, if applicable.
2. Lift Algidex Ag I.V. Patch from the catheter and discard according to institutional policy.
3. Gently wipe any Algidex matrix residue during routine dressing change using saline or antiseptic skin prep as stated in the institutional procedural guidelines for vascular dressing changes.

Product #	Pad	Qty/Bx	Qty/Cs
46-IV22	1" disc	5	50
46-IV25	1" disc w/7mm opening	5	50
46-IV375	1½" disc w/7mm opening	5	50



200 DeBusk Lane  
Powell, Tennessee 37849  
[www.deroyal.com](http://www.deroyal.com)

Reprint #0-2043  
Issue date: 3/06  
©2005 DeRoyal