EN: DEROYAL® RIGID LIGHT HANDLE

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NL: DEROYAL® RIGIDE LAMPHANDVAT

MANUFACTURED BY:
DeRoyal Industries, Inc.
200 DeBusk Lane
Powell, TN 37849 USA
888.938.7828 or (001) 865.938.7828
www.deroyal.com

MedEnvoy Global B.V.
Prinses Margrietplantsoen 33 – Suite 123
2595 AM the Hague
The Netherlands

PART#4-10702G | REVISED 6/24

©2024 DeRoyal Industries, Inc.

All Rights Reserved. DeRoyal and the DeRoyal logo are trademarks or registered trademarks of DeRoyal Industries, Inc.

DEROYAL® RIGID LIGHT HANDLE

	DO NOT REUSE
	MEDICAL DEVICE
	DO NOT USE IF PACKAGE IS DAMAGED
	NOT MADE WITH NATURAL RUBBER LATEX
	FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.
STERILE PRODUCTS ONLY	
	STERILIZED WITH ETHYLENE OXIDE
	DO NOT RESTERILIZE
	SINGLE STERILE BARRIER SYSTEM

IMPORTANT INFORMATION
Please read all instructions, warnings, and cautions before use. Correct application is essential for proper functioning of product.

INTENDED USE
The DeRoyal® Rigid Light Handle Cover is intended to screw-in to an adapter for surgical lights that allows for light manipulation and adjustment while maintaining a sterile field.

INTENDED USERS
The DeRoyal® Screw-in Rigid Light Handle Cover is intended to be used only by a licensed practitioner or trained healthcare professionals.

PATIENT TARGET GROUP(S)
Device is non-patient contacting.

INDICATIONS
Use this product whenever surgical lights are required during a surgical case and will need to be adjusted by personnel who are operating within the sterile field.

WARNINGS

- This product should only be used on surgical lights where an appropriate adapter has been securely installed. The adapter should be permanently installed and should not be removed between surgical cases.
- DO NOT** reuse, reprocess, or re-sterilize. This device is for single use only. Reuse, reprocessing, or re-sterilization may damage the device, and may result in user or patient injury, patient infection, or patient death.
- Product should not be used if sterile packaging is damaged or opened. If damage is found, remove the product from the operative area to prevent unintended use.
- DO NOT** use if product appears visually damaged or degraded.

PRECAUTIONS

- The handle should be screwed in tightly and securely into the surgical lights by hand only.
- DO NOT** over-tighten the handle.

- DO NOT** use pliers, wrenches, or other tools to tighten the handle.
- Verify securement of the light handle to the adaptor.

CAUTIONS

- Care should be used in not hitting one's head on the handle. The handle will protrude at least 4 inches (10.2 cm) beyond the center of the surgical light. This could cause personal injury and will contaminate the sterile field.

DIRECTIONS FOR USE

- Visually examine the package for damage before opening the package. Replace the light handle(s) if the package or the handle(s) themselves show signs of damage.
- Open package using aseptic technique.
- Separate the light handles from each other by sliding them parallel to each other.
- Squeeze the flange together on either side of the shaft in order to release the flange tabs from the shaft of the handle(s).
- Invert the flange such that it creates a concave "dish" around the threaded portion of the handle(s).
- Screw the handle(s) clockwise into the surgical light handle socket or the adapter permanently attached to the light. Screw the handle(s) in hand-tight.
- If center focus lights are used, be sure to screw the handle in clockwise the full rotation of the center focus adjustment.
- After use, unscrew the handle(s) counterclockwise and remove from the surgical light.
- Discard the handle(s) after use. Do not re-sterilize.

STERIZATION INFORMATION
NON-STERILE PRODUCT: DeRoyal intends that non-sterile product for use in sterile environments be further processed by our customers, including further packaging and/or sterilization according to their own validated processes. These products are to be sterilized by ethylene oxide only, as other sterilization methods have not been validated by DeRoyal and may damage the product which could cause a device malfunction and/or injury to the patient.

DISPOSAL
Dispose of any used device according to local, state, and federal laws and regulations. For safe disposal of devices, follow your facility's protocol.

	KEEP DRY
	KEEP AWAY FROM SUNLIGHT

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

WARRANTY
DeRoyal products are warranted for one hundred and twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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