EN: FOR USE WITH YSI 400 SERIES ELECTRONICS ONLY
ES: PARA USO SOLO CON COMPONENTES ELECTRÓNICOS SERIE YSI 400
FR: À N'UTILISER QU'AVEC LE MATÉRIEL ÉLECTRIQUE SÉRIE YSI 400
PT: PARA USO SOMENTE COM APARELHOS ELETRÔNICOS SÉRIE YSI 400

DeRoyal uses the symbols and meanings from standards ISO 15523-1. These symbols are placed next to the text explaining their meaning in this instructions for use (IFU), a complete symbols glossary is available online at deroyal.com/symbols or by contacting customer service.

Manufactured by:
 DeRoyal Industries, Inc.
 200 DeBusk Lane
 Powell, TN 37849 USA
 888.938.7828 or (001) 865.938.7828
 www.deroyal.com

PART#74-14905
 Revised 2/2025
 ©2025 DeRoyal Industries, Inc.
 All Rights Reserved. DeRoyal and the DeRoyal logo are registered trademarks of DeRoyal Industries, Inc.

ESOPHAGEAL STETHOSCOPE WITH TEMPERATURE SENSOR

	MR UNSAFE
	DO NOT RESTERILIZE
	DO NOT REUSE
	STERILIZED WITH ETHYLENE OXIDE
	NOT MADE WITH NATURAL RUBBER LATEX
RX ONLY	FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.

INDICATIONS

The DeRoyal® Esophageal Stethoscope with Temperature Sensor is indicated to be used for routine monitoring of core body temperature, as well as to monitor heart and respiratory sounds in an anesthetized patient by inserting the stethoscope tube into the esophagus.

CONTRAINDICATIONS

The use of an esophageal stethoscope may be contraindicated in patients with esophageal diverticulum or stenosis, in neonates or small infants undergoing neck surgery, in patients undergoing a tracheostomy or insertion of an internal jugular catheter, and in patients with conditions deemed unsafe for insertion of the stethoscope into the esophagus.

! WARNINGS

- Tracheal or bronchial insertion of stethoscope may cause airway obstruction and may be fatal if uncorrected. Correct placement may be verified with a laryngoscope.
 - Carefully inspect the packaging and device for damage or defect. **DO NOT** use if the packaging is open or the device is damaged.
- This is a single use device. It is not intended for disinfection and/or subsequent re-use. Re-use could result in microbial contamination and potential patient injury.
- This device has not been evaluated for reprocessing or re-sterilization. Reprocessing and/or re-sterilization may damage the device, rendering it unusable and/or may lead to device failure.

ADVERSE REACTIONS

While adverse reactions are rare, the following are potential adverse reactions associated with the use of temperature probes during insertion or while the device is in use:

- infection/sepsis
- allergic reaction to materials
- airway obstruction
- aspiration pneumonitis
- bronchial insertion
- electrical burns and shock
- epistaxis
- esophageal abrasion
- esophageal perforation
- tracheal insertion
- trauma to pharynx

CAUTIONS

- Use of improper insertion or removal technique may damage the device causing foreign body retention.
- Improper use or positioning of the device can produce inaccurate temperature readings that could result in

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patient hyperthermia, hypothermia, or an inaccurate diagnosis.

- If incorrect temperature readings or excessive response time is noticed, discard probe and replace.
- Intubate patient with tracheal tube before inserting the esophageal stethoscope.
- Inform the surgeon that an esophageal stethoscope has been inserted.
- Lubricate the stethoscope before insertion and use accepted medical techniques in the insertion and removal of the stethoscope.
- Take care to ensure that the cable and connector do not get wet and proper techniques are used during electrosurgical procedures to reduce radio frequency interference current and potential burn to the patient.
- An adequate electrosurgical dispersive ground electrode close to the active surgical site should be properly connected.
- Do not intertwine the cables, especially the monitor cables with the electrosurgical unit's cables.
- The operation of the patient temperature monitor may be temporarily affected during electrosurgical activations. Unusual temperature readings should be checked.
- If unusual or inconsistent temperature readings are noticed, check connections. If issue is unresolved, replace temperature probes.
- Check the temperature monitoring equipment user manual for information regarding compliance with Clause 5 of IEC 60601-1-2 - Medical electrical equipment – Part 1-2. General requirements for basic safety and essential performance – Collateral Standard: Electromagnetic disturbances – Requirements and tests.
- The probe is intended to be used with a DeRoyal-branded interface cable. Use with an incompatible cable may affect performance.

DIRECTIONS FOR USE

1. Remove the Esophageal Stethoscope with Temperature Sensor from its sterile package using aseptic technique.
2. Verify the compatibility of the probe, interface cable, and patient monitor before use.
3. Lubricate the cuff with a suitable water soluble lubricant.
4. Insert the stethoscope and guide it into the esophagus.
5. Listen for airway leaks when inserting. Invasion of the trachea may permit air to escape around the tracheal tube cuff. Do not seal the leak. Remove the esophageal stethoscope and reinsert. To verify proper placement, use a laryngoscope.
6. Adjust the stethoscope's depth in the esophagus to maximize the heart and respiratory sounds.
7. Connect the temperature probe to the appropriate reusable temperature cable from DeRoyal®. Connect the cable to the patient monitor. Secure the cable using the drape clip.

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EN: ESOPHAGEAL STETHOSCOPE WITH TEMPERATURE SENSOR

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8. Follow the directions for use of the patient temperature monitor instrument.
9. Remove the probe prior to extubation of the tracheal tube following accepted medical technique.
10. Disconnect the probe (at the connector) and discard the probe. To disconnect, grasp both connectors firmly and pull. **DO NOT** pull cable or wire.

! SAFETY IN MRI NOT EVALUATED

The Esophageal Stethoscope with Temperature Sensor's safety in MRI has not been evaluated, hence, the probe is considered MR Unsafe.

OPERATING SPECIFICATIONS

- Operating Temperature: 25° C to 45° C
- Rated Output Range: 35° C to 42° C
- Rated Extended Output Range: 25° C to 35° C and 42° C to 45° C
- The reference body site is the core body temperature.
- The measuring site is the esophagus.
- Accuracy in Rated Output Range: +/- 0.2° C
- Accuracy in Rated Extended Output Range: +/- 0.3° C
- Time response
 - The temperature probe has a direct mode of operation.
 - The maximum heating transient time is 2 minutes 45 seconds.
 - The maximum cooling transient time is 1 minute 45 seconds.

DISPOSAL

Dispose of any used device according to local, state, and federal laws and regulations. For safe disposal of devices, follow your facility's protocol.

STORAGE AND TRANSPORT CONDITIONS

	KEEP DRY
	KEEP AWAY FROM SUNLIGHT

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

WARRANTY

DeRoyal® products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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