GENERAL PURPOSE TEMPERATURE PROBE

SONDA DE TEMPERATURA MULTIUSO

SONDE DE TEMPÉRATURE UNIVERSELLE

SONDA DE TEMPERATURA PARA FINALIDADES GERAIS

Manufactured by:
DeRoyal Industries, Inc.
200 DeBusk Lane
Powell, TN 37849 USA
888.938.7828 or (001) 865.938.7828
www.deroyal.com

Part#74-14906

Revised 2/2021

©2021 DeRoyal Industries, Inc.

All Rights Reserved. DeRoyal, the DeRoyal logo and Improving Care. Improving Business. is a registered trademark of DeRoyal Industries, Inc.

Mindray is a registered trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Draeger is a registered trademark of Draegerwerk AG & CO. KGaA Corporation.

GENERAL PURPOSE TEMPERATURE PROBE

	DO NOT RESTERILIZE
	DO NOT REUSE
	STERILIZED WITH ETHYLENE OXIDE
	NOT MADE WITH NATURAL RUBBER LATEX
	FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.

IMPORTANT

Please read all warnings and instructions before use. Correct application is essential for proper product function and to reduce the risk of injury.

INTENDED USE

The DeRoyal® General Purpose Temperature Probe is to be used for routine monitoring of core body temperature.

PRODUCT DESCRIPTION

The General Purpose Temperature Probe consists of a temperature-sensing wire set that is enclosed within a small-diameter plastic tube for ease of application. The probe is available in 9 and 12 French sizes and may be inserted into the nasopharyngeal, esophageal, and rectal cavities.

The temperature monitoring sensor is located at the tip of the probe and secured inside the lumen of the tubing. The blue connector at the proximal end interfaces with DeRoyal-approved cables for connection to a monitoring unit to display temperature. The temperature monitoring sensor is available in YSI 400 or YSI 700 series configurations for connection to a compatible monitoring unit. Whether the probe is a 400 or 700 series configuration is indicated on the device labeling. Consult the user manual for the patient monitor to ensure compatibility. The probes are designed to interface with DeRoyal-branded cables for connection with YSI 400 or 700 series compatible monitors, including the following patient monitors and equivalent models: Mindray® Passport, Philips IntelliVue®, Siemens/Draeger Infinity®, and GE Datex-Ohmeda® brands.

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CONTRAINDICACIONES

CONTRAINDICACIONES

The use of an esophageal probe may be contraindicated in patients with esophageal diverticulum or stenosis, in neonates or small infants undergoing neck surgery, and in patients undergoing a tracheostomy or insertion of an internal jugular catheter. Rectal placement is contraindicated in neonates and small infants due to possibility of colon perforation. Nasopharyngeal insertions not recommended for patients who are taking anticoagulants or when manipulation of mucosa is undesirable.

! WARNINGS

- Tracheal or bronchial insertion of the probe may cause airway obstruction and may be fatal if uncorrected. Correct placement may be verified with a laryngoscope.
- While adverse reactions are rare, the following have been reported to be associated with the use of temperature probes during insertion or while the device was in use:
 - airway obstruction
 - aspiration pneumonia
 - bronchial insertion
 - electrical burns
 - epistaxis
 - esophageal abrasion
 - esophageal perforation/colon perforation
 - rectal bleeding
 - stool insertion
 - tracheal insertion
 - trauma to pharynx
- May cause epistaxis (nosebleed). In the pregnant patient, children with large adenoids, and patients with clotting disorders, the incidence of epistaxis is high.
- The nasal mucosal temperature can vary based on exact location (i.e., anterior versus posterior mucosa) and respiratory cycle (i.e., inspiration versus expiration). The exact location and timing of the respiratory cycle should be noted and taken into consideration when using a nasopharyngeal or nasal temperature probe.
- Inspect packaging for any damage and product for any damage, contaminate or missing parts prior to use.

CAUTIONS

- If the probe is to be used in the esophagus on an

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anesthetized patient, intubate the patient with the endotracheal tube before inserting the probe.

- Lubricate the probe with a medical grade water-soluble lubricant before insertion and use accepted medical techniques during insertion and removal of the probe.
- Take care to ensure that the cable and connector do not get wet and that proper techniques are used during electrosurgical procedures to reduce radio frequency interference current and potential burns to the patient.
- An adequate electrosurgical dispersive ground electrode close to the active surgical site should be properly connected.
- DO NOT** intertwine the cables, especially the monitor cables, with the electrosurgical unit's cables.
- The operation of the patient temperature monitor may be temporarily affected during electrosurgical activations. Unusual temperature readings should be checked.
- If probe has inaccurate, unstable, or no temperature readings, discard and replace.
- The probe is intended to be used with DeRoyal-approved interface cables. Use with an incompatible cable may affect performance.

MR COMPATIBILITY

- The DeRoyal® General Purpose Temperature Probe has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the GP Probe in the MR environment is unknown. Scanning a patient who has this medical device may result in patient injury.

DIRECTIONS FOR USE

- Remove the probe from its sterile package.
- Verify the compatibility of the probe, interface cable, and patient monitor.
- Lubricate the probe tip with a suitable water soluble lubricant.
- Insert the probe and guide it into either the rectum, esophagus, nostrils or nasopharynx cavity.
- Connect the temperature probe to the appropriate reusable temperature cable from DeRoyal. Connect the cable to the patient monitor. Secure the cable using drape clip.
- Follow the directions for use of the patient

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temperature monitor instrument.

- If used in the esophagus, remove the probe prior to extubation of the tracheal tube following accepted medical technique.
- Disconnect the probe (at the connector) and discard the probe. To disconnect, grasp both connectors firmly and pull. **DO NOT** pull on cable or wire.
- This probe is **NOT** reusable.
- Dispose of probe according to facility protocol.

OPERATING SPECIFICATIONS

- Rated output range: 25° to 45° Celsius
- The reference body site is the core body temperature.
- The measuring site(s) are Esophagus, Rectum, or Nasopharynx.
- Accurate to +/- 0.2° Celsius
- Time response:
 - The temperature probe has a direct mode of operation.
 - The maximum heating transient time is 1 minute
 - The maximum cooling transient time is 30 seconds

STORAGE AND TRANSPORT CONDITIONS

	KEEP DRY
	KEEP AWAY FROM SUNLIGHT
	DO NOT USE IF PACKAGE IS DAMAGED
	RELATIVE HUMIDITY: 85% NON-CONDENSING
	STORAGE TEMPERATURE: - 25°C TO + 55°C
	OPERATING TEMPERATURE: 25°C TO 45°C

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

WARRANTY

DeRoyal® products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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