another steel clamp and incise the cord between the clamps.

- Using the steel clamp, lower the section of cord into the top of the Umbilicup™ device.
- Slowly release the steel clamp, allowing umbilical blood to accumulate in the collection funnel. Place cap on the top of the device.
- Insert empty vacuum collection tube into needle sleeve at bottom of the Umbilicup™ device.
- Once filled, remove collection tube from the device and transport to lab for designated tests.
- Discard into the facility's sharps disposal container.

#### DIRECTIONS FOR CORD BLOOD GAS TRANSPORT

- Notify the laboratory that blood gas specimen is being collected.
- After delivery, double-clamp the cord; use a plastic clamp above the Pean or Kocher forceps. Cut the cord.
- Place a second plastic clamp 2 to 3.7 inches beyond the first clamp and a stainless steel clamp just above the second clamp. Cut the cord again between these clamps.
- The cord segment, secured by two plastic clamps, is placed into the top of the Umbilicup™ device.
- Pack sterile ice into the top of the device. Place the cap on the top of the Umbilicup™ device.
- Promptly transport the Umbilicup™ device containing the section of umbilical cord to the laboratory.

#### STERILIZATION INFORMATION

DeRoyal intends that non-sterile products for use in sterile environments be further processed by our customers, including further packaging and/or sterilization. These products are to be sterilized by ethylene oxide, utilizing the following parameters: Temperature: 120°F - 135°F (48.9°C - 57.2°C), Ethylene Oxide Gas Concentration: 625 - 800 mg/l, Exposure Time: 4 hrs minimum (Preconditioning: Temperature 100°F -120°F (37.8°C - 48.9°C), Humidity 55-75%, Time 18hrs minimum), (Aeration: Temperature 100°F -135°F (37.8°C - 42.2°C), Time 18hrs minimum). Other sterilization methods have not been validated by DeRoyal and may damage the product which could cause a device malfunction and/or injury to the patient.

The device has not been evaluated for re-processing or re-sterilization. Reprocessing and/or re-sterilization may damage the device, rendering it unusable and/or may lead to device failure which could result in patient illness, injury or death.

#### STORAGE AND TRANSPORT CONDITIONS

Do not store product at extreme temperatures or in a moist/damp environment.

	KEEP DRY
	KEEP AWAY FROM SUNLIGHT

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

#### WARRANTY

DEROYAL® PRODUCTS ARE WARRANTED FOR NINETY (90) DAYS FROM THE DATE OF SHIPMENT FROM DEROYAL AS TO PRODUCT QUALITY AND WORKMANSHIP. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

Improving Care. Improving Business.®



Manufactured by:  
DeRoyal Industries, Inc.  
200 DeBusk Lane  
Powell, TN 37849 USA  
888.938.7828 or (001) 865.938.7828  
www.deroyal.com

Part#0-10569 Revised 8/2022

©2022 DeRoyal Industries, Inc.  
All Rights Reserved. DeRoyal, the DeRoyal logo, Improving Care. Improving Business., and Umbilicup are trademarks or registered trademarks of DeRoyal Industries, Inc.

#### UMBILICUP™ UMBILICAL CORD BLOOD COLLECTION DEVICE

	DO NOT REUSE
	STERILIZED WITH ETHYLENE OXIDE
	NOT MADE WITH NATURAL RUBBER LATEX
	FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.

#### IMPORTANT INFORMATION

 Please read all instructions, warnings, and precautions before use. Correct application is essential for proper functioning of product.

#### INTENDED USE

The DeRoyal® Umbilicup™ Umbilical Cord Blood Collection Device is intended to be used by a licensed practitioner for the collection and transport of umbilical cord blood for the purposes of routine laboratory testing.

#### INDICATIONS


Use of this device is indicated when it is necessary to collect umbilical cord blood for laboratory testing of newborns after either vaginal or cesarean section deliveries. The laboratory tests on the umbilical cord blood may be performed depending on local law or on the health and condition of the mother or baby as deemed and hereto ordered by the physician.

#### CONTRAINDICATIONS

The Umbilicup™ device is not intended for the collection, transportation, or storage of umbilical cord blood intended for stem cell research or banking.

#### ⚠ CAUTIONS

- Product should not be used if the cup is broken, cracked or chipped. Product should not be used if the cap or lid does not fit. If damage is found, remove the product from the sterile procedure area to prevent unintended use and return to the manufacturer.
- The package should be opened using sterile technique.
- Product should not be used if sterile packaging is damaged or opened or if shelf-life has expired. If damage is found or the product is expired, remove the product from the sterile procedure area to prevent unintended use and return to the manufacturer.
- This device is for single use only. It is not intended for disinfection and/or subsequent re-use, which can result in microbial contamination causing health deterioration of the patient or inaccurate laboratory results.

- DO NOT** re-sterilize. This device has not been evaluated for reprocessing or re-sterilization. Reprocessing and/or re-sterilization may damage the device, rendering it unusable and/or may lead to device failure, which could result in patient illness, injury or death.

#### WARNINGS

- Care should be exercised whenever handling blood to prevent splashes, sprays, and spills. Proper personal protective equipment should be worn according to local or facility policies.
- Care should be exercised whenever handling sharp needles to prevent accidental needle sticks.
- Product should not be used if the needle is detached, loose, or exposed.
- If needle becomes detached with the vacuum tube is inserted or removed, or if the needle sheath does not return to cover the needle, caution should be exercised when handling the device or separated needle.
- Product should be disposed of in a facility approved sharps container.

#### DIRECTIONS FOR GENERAL BLOOD COLLECTION METHOD #1 - GRAVITY FLOW METHOD

- Double clamp and cut umbilical cord as usual.
- Place clamped free-end of the cord into Umbilicup™ device.
- Release cord clamp slowly to allow gravity flow to drain cord blood into the device.
- Insert empty vacuum collection tube into sheath-protected needle at bottom of the Umbilicup™ device.
- Once filled, remove collection tube from the device and send to the lab for designated studies.
- Place lid firmly on top of the Umbilicup™ device and place red cap onto bottom of device by pressing down on a flat surface.
- Discard into the facility's sharps disposal container.

#### DIRECTIONS FOR GENERAL BLOOD COLLECTION METHOD #2 - DOUBLE CLAMP METHOD

- Double clamp and cut umbilical cord as usual.
- Place plastic cord clamp approximately 6-7 inches (15-17 cm) above the metal clamp. Immediately above the plastic clamp, apply

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