



# WARRIOR<sup>®</sup> RECOVERY

Post-Operative  
Knee Brace



# WARRIOR<sup>®</sup>

## RECOVERY

### Post-Operative Knee Brace

- **Easy-to-use**, lever-action, micro-adjustable telescoping struts for length adjustments from 18" to 29"
- **Four convenient preset** drop lock settings, 0°, 15°, 30°, 45°
- **Enhanced range-of-motion (ROM)** hinge interface for easy setting and patient compliance
- **ROM** is set in 10° increments
- **Soft, circumferential and trimmable straps** plus removable Y-tab feature optimizes individual patient fit and comfort
- **Inner pad panel adjustment** (2mm hex wrench, not included) to avoid incision sites and provide a more customized fit
- **Customizable anti-migration strap pads** for superior fixture and reduced slippage
- **Integrated "bend zone"** on lower hinge arms allows for customization and accommodation to varying patient anatomy
- **Laminate reinforced padding** to support prolonged use and maintain pad integrity
- **Extended padding** at malleoli to prevent irritation and improve comfort during use

#### Indications:

Post-operative stabilization and rehabilitation; Post total knee replacement; Post ORIF; Fractures of the femur, tibia, fibula or patella; Ligament reconstructive procedures; Meniscal repairs; Post-fracture management.

Size	Description	Adjustment Length	Product Number
Universal	Warrior <sup>®</sup> Recovery Post-Op Knee Brace	18" to 29"	KB9000-01

Warrior<sup>®</sup> is a registered trademark of DeRoyal Industries, Inc.



Improving Care. Improving Business.<sup>®</sup>

