



Improving Care. Improving Business.®

PRO-II® NPWT DEVICE

EXTERNAL EXHAUST FILTER CHANGE SET

NON-STERILE

For use only with PRO-II® Negative Pressure Wound Therapy devices.

PRODUCT DESCRIPTION:

The external exhaust filter change set for the Prospera® PRO-II device contains parts and tools needed to properly change the exhaust filter on the device.

Ensure all parts are included in the kit prior to starting the replacement.

PARTS LIST:

- 1 Charcoal Exhaust filter
- 2 PRO-II housing feet (one needed for replacement, 1 extra provided)
- 1 2MM Allen Wrench
- 1 Disinfecting Wipe

You will need a flat head screw driver or needle nose pliers, not provided.

REMOVING THE USED FILTER

1. Ensure the PRO-II device is not connected to a power source, is completely shut off and no canister is attached to the device.
2. Flip the device over so the top of the device (screen side) is face down on a flat surface.
3. Remove the clear rubber foot BELOW the outlet air (ABLUF) arrow **only** using a flat head screwdriver or needle nose pliers.
4. Once removed, use the allen wrench provided to remove the screw by rotating counterclockwise.
5. The lid over the exhaust compartment should come out easily once the screw is removed. PLACE THE SCREW IN A SECURE LOCATION.
6. Remove the filter. Use allen wrench or screwdriver to loosen if necessary.

INSERTING NEW EXHAUST FILTER

1. Once old filter is removed, if necessary, wipe with disinfecting wipe provided.
2. Allow filter compartment to dry fully prior to inserting new charcoal filter.
3. Insert new charcoal filter. It is dye-cut to fit exactly. **Note:** there is no top/bottom.
4. Place the plastic cover back on the bottom of the

- housing.
- 5. Replace original screw and fasten using the provided allen wrench by rotating clockwise until snug. **DO NOT OVERTIGHTEN.**
- 6. Replace clear rubber foot with new foot provided.
- 7. Perform function test as per facility protocols.

WARRANTY

DEROYAL PRODUCTS ARE WARRANTED FOR NINETY (90) DAYS FROM THE DATE OF SHIPMENT FROM DEROYAL AS TO PRODUCT QUALITY AND WORKMANSHIP. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**



Improving Care. Improving Business.®

Manufactured for DeRoyal Industries, Inc.
200 DeBusk Lane | Powell, TN 37849 U.S.A.
888.938.7828 | 865.938.7828 | www.deroyal.com

Part#1003340 Revision 9/16
©2016 DeRoyal Industries, Inc. All Rights Reserved.
DeRoyal, the DeRoyal logo, PRO-II, Prospera, and Improving Care. Improving Business. are trademarks or registered trademarks of DeRoyal Industries, Inc.



Improving Care. Improving Business.®

PRO-II® NPWT DEVICE

EXTERNAL EXHAUST FILTER CHANGE SET

NON-STERILE

For use only with PRO-II® Negative Pressure Wound Therapy devices.

PRODUCT DESCRIPTION:

The external exhaust filter change set for the Prospera® PRO-II device contains parts and tools needed to properly change the exhaust filter on the device.

Ensure all parts are included in the kit prior to starting the replacement.

PARTS LIST:

- 1 Charcoal Exhaust filter
- 2 PRO-II housing feet (one needed for replacement, 1 extra provided)
- 1 2MM Allen Wrench
- 1 Disinfecting Wipe

You will need a flat head screw driver or needle nose pliers, not provided.

REMOVING THE USED FILTER

1. Ensure the PRO-II device is not connected to a power source, is completely shut off and no canister is attached to the device.
2. Flip the device over so the top of the device (screen side) is face down on a flat surface.
3. Remove the clear rubber foot BELOW the outlet air (ABLUF) arrow **only** using a flat head screwdriver or needle nose pliers.
4. Once removed, use the allen wrench provided to remove the screw by rotating counterclockwise.
5. The lid over the exhaust compartment should come out easily once the screw is removed. PLACE THE SCREW IN A SECURE LOCATION.
6. Remove the filter. Use allen wrench or screwdriver to loosen if necessary.

INSERTING NEW EXHAUST FILTER

1. Once old filter is removed, if necessary, wipe with disinfecting wipe provided.
2. Allow filter compartment to dry fully prior to inserting new charcoal filter.
3. Insert new charcoal filter. It is dye-cut to fit exactly. **Note:** there is no top/bottom.
4. Place the plastic cover back on the bottom of the

- housing.
- 5. Replace original screw and fasten using the provided allen wrench by rotating clockwise until snug. **DO NOT OVERTIGHTEN.**
- 6. Replace clear rubber foot with new foot provided.
- 7. Perform function test as per facility protocols.

WARRANTY

DEROYAL PRODUCTS ARE WARRANTED FOR NINETY (90) DAYS FROM THE DATE OF SHIPMENT FROM DEROYAL AS TO PRODUCT QUALITY AND WORKMANSHIP. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**



Improving Care. Improving Business.®

Manufactured for DeRoyal Industries, Inc.
200 DeBusk Lane | Powell, TN 37849 U.S.A.
888.938.7828 | 865.938.7828 | www.deroyal.com

Part#1003340 Revision 9/16
©2016 DeRoyal Industries, Inc. All Rights Reserved.
DeRoyal, the DeRoyal logo, PRO-II, Prospera, and Improving Care. Improving Business. are trademarks or registered trademarks of DeRoyal Industries, Inc.