




Oil Emulsion



Xeroform

Xeroform and Oil Emulsion are sterile dressings used to treat everything from minimal drainage to moderate and heavy drainage wounds such as surgical incisions, lacerations and burns to skin grafts, donor sites, dermal ulcers, and more. They are non-adherent and help maintain a moist wound environment.



## OIL EMULSION

non-adhering wound dressing



- A sterile, knitted cellulose fabric dressing impregnated with oil emulsion blend
- Non-adherent to wound site
- Non-occlusive, allowing flow of exudate away from wound site into secondary dressing
- Helps maintain a moist wound environment
- Can be cut to shape or size



## XEROFORM

petrolatum wound dressing



- A sterile, fine mesh gauze impregnated with a blend of 3% Bismuth Tribromophenate (Xeroform) and USP Petrolatum
- Non-adherent to wound site
- Helps maintain a moist wound environment
- 3% Bismuth Tribromophenate (Xeroform) provides deodorizing action

Product #	Size	Qty/Bx	Bx/Cs	Qty/Cs	HCPC
46-620	3" x 3"	50	12	600	A6222
46-621	3" x 8"	24	6	144	A6223
46-622	3" x 16"	36	6	216	A6223


Product #	Size	Qty/Bx	Bx/Cs	Qty/Cs	HCPC
46-610	1" x 8"	50	4	200	A6222
46-611	5" x 9"	50	4	200	A6223
46-612	4" x 3 yds	6	6	36	A6266

### indications:

Surgical incisions  
Burns  
Skin grafts  
Donor sites  
Dermal ulcers  
Suture lines  
Lacerations

### indications:

Minimally draining wounds  
Surgical incisions  
Lacerations  
Donor sites  
Burns



\* Any L-Code or ICD-10 Code information provided by DeRoyal is only a recommendation provided for the convenience of and as a resource to healthcare professionals. These recommendations are based upon private and public information including L-Code descriptors, ICD-10 diagnostic codes, Medicare PDAC information and information from experienced consultants and DMEPOS suppliers. DeRoyal does not guarantee or imply claim reimbursement or coverage for any specific DeRoyal product or in any specific amount. For more information, please visit our website.

DeRoyal and DeRoyal logo, Improving Care. Improving Business. are registered trademarks of DeRoyal Industries, Inc.

© 2016, DeRoyal. All rights reserved.

DeRoyal Industries, Inc. | 200 DeBusk Lane, Powell, TN 37849 USA

800.DEROYAL or 865.938.7828 | www.deroyal.com

Reprint # 0-1885 | Rev. 11/16

DeRoyal®  
Improving Care. Improving Business.™