




Improving Care. Improving Business.



DeRoyal presents the revolutionary  
**Element<sup>®</sup> Sport**

An ankle brace engineered to provide true subtalar control through calcaneal stabilization and positioning. The low profile design fits comfortably in athletic footwear, helping to increase patient compliance. The unique patented and patent pending heel control strapping system holds the calcaneus under the talus; thus, controlling subtalar inversion and talar rotation to assist in tissue healing, preventing chronic instability and allowing for a quicker return to activity. In addition, the **Element<sup>®</sup> Sport** ankle brace provides standard stirrup stabilization to reduce ankle joint inversion and eversion.



THE NEW GENERATION OF CONTROL, COMFORT & SUPPORT

888.938.7828 | 865.938.7828 | WWW.DEROYAL.COM

REV 10/10 REPRINT# 0-2214