

WARRIOR®



COMFORT for the way you **LIVE**

CONTROL for the way you **PLAY**

CONFIDENCE in how you **RECOVER**



Improving Care. Improving Business.®



FEATURES & BENEFITS

- Classic pull-up design, wrap-around design, and new half-wrap design to optimize individual patient fit and comfort
- Integrated anti-migration strip (pull-up designs) for superior fixture and reduced slippage
- Flex-tab feature to prevent hinge digging or irritation
- Plaited stretch-nylon material (engineered for swimwear) helps wick away moisture and promote aeration behind the leg and popliteal region
- Waffled spacer fabric promotes breathability and a cooler environment
- Removable condyle pads for increased comfort and stability
- Easy-grip silicone pull tabs for easier donning and doffing
- Two detachable, medial/lateral half-moon shaped buttresses (pull-up style) and one detachable, inferior half-moon shaped buttress (half-wrap) for patellar control
- Optional flexion/extension stops (0° to 100° in 10° increments) for range motion control
- Optional “hyper control” strapping to help prevent hyperextension (regular designs)



INDICATIONS: Mild to moderate ACL, PCL, MCL, and LCL sprains
 Combined ligament instabilities | Patella subluxation

SIZE	MEASUREMENT*	PULL-UP LONG	PULL-UP SHORT	WRAP AROUND LONG	WRAP AROUND SHORT	HALF WRAP LONG	HALF WRAP SHORT
X-SMALL	13" - 15½"	14011004	14021004	14012004	14022004	14013004	14023004
SMALL	15½" - 18½"	14011005	14021005	14012005	14022005	14013005	14023005
MEDIUM	18½" - 21"	14011006	14021006	14012006	14022006	14013006	14023006
LARGE	21 - 23½"	14011007	14021007	14012007	14022007	14013007	14023007
X-LARGE	23½" - 26½"	14011008	14021008	14012008	14022008	14013008	14023008
XX-LARGE	26½" - 29½"	14011009	14021009	14012009	14022009	14013009	14023009
XXX-LARGE	29½" - 32"	14011010	14021010	14012010	14022010	14013010	14023010

*Measurement should be taken 6" above mid-patella



Improving Care. Improving Business.®

DeRoyal Industries, Inc.
 200 DeBusk Ln | Powell, TN 37849
 888.938.7828 | www.deroyal.com
 Revision 2/14 | Reprint # 0-2339

