

Capital™ Cervical Collar



 DeRoyal®

Improving Care. Improving Business.®

Capital™ Cervical Collar



FEATURES & BENEFITS

- Designed to help restrict flexion, extension, lateral bending, and axial rotation
- Reinforced sternal padding for optimum support and comfort
- Shoulder, mandible, and occipital contours for increased stability
- Safe for use with X-Ray, CT, and MRI
- Sheer finished open-cell liner is specifically designed to help prevent skin breakdown and provide sustained comfort during long-term use
- Replacement pad set is available upon request
- Not made with natural rubber latex

INDICATIONS

Cervical sprains and strains | Cervical spinal stenosis | Cervical disc displacement
 Cervical arthritis | Whiplash injury | Cervical spondylosis and spondylolisthesis
 Degenerative disc disease | Cervicalgia

| Description | Size | Collar | Collar w/Extra Pad Set | Pad Set |
|-------------------------|--------------------|----------|------------------------|----------|
| Capital Cervical Collar | Super Short | 11421000 | 11422000 | 11423000 |
| Capital Cervical Collar | Stout | 11431000 | 11432000 | 11433000 |
| Capital Cervical Collar | XS | 11441000 | 11442000 | 11443000 |
| Capital Cervical Collar | Short | 11451000 | 11452000 | 11423000 |
| Capital Cervical Collar | Regular | 11461000 | 11462000 | 11423000 |
| Capital Cervical Collar | Large | 11471000 | 11472000 | 11423000 |
| Capital Cervical Collar | 0-6 months | 11481000 | 11482000 | 11483000 |
| Capital Cervical Collar | 6 months - 2 years | 11481001 | 11482001 | 11483001 |
| Capital Cervical Collar | 2-6 years | 11481002 | 11482002 | 11483002 |
| Capital Cervical Collar | 6-12 years | 11481003 | 11482003 | 11483003 |



Improving Care. Improving Business.®

DeRoyal Industries, Inc.
 200 DeBusk Lane | Powell, TN 37849
 888.938.7828 or 865.938.7828
www.deroyal.com

