



Your Frontline Defense for Optimum Fluid Control



Improving Care. Improving Business.®

Premium Solidifier+™

fluid solidifier plus, treatment technology



- Not only solidifies liquid medical waste, but is also a U.S. Environmental Protection Agency (EPA) registered liquid medical waste treatment
- Treats the contents of suction canisters and liners rendering them non-infectious and suitable for white bag disposal in most states¹
- Works like most standard solidifiers; except that when the suction canister is solidified, it is also treated and ready for disposal

Part Number	Description	Case Quantity
71-1000S	1000cc Treatment Technology	100 Per Case
71-1500S	1500cc Treatment Technology	100 Per Case
71-3000S	3000cc Treatment Technology	100 Per Case


Premium Solidifier Snapshot:

- Quickly encapsulates large volumes of bodily fluids
- One ounce cost-effectively solidifies 1000cc of liquid
- Works well with all suction canisters and liners
- Eliminates transport risk of biohazardous fluids
- Eliminates time and risk when emptying the collection unit
- Not made with natural rubber latex

Premium Solidifier+™

fluid solidifier plus, treatment technology

States approved for White Bag Disposal



¹ **White Bag Approved** - Alabama, Alaska, Arizona, Arkansas, Colorado, Connecticut, Delaware, District of Columbia, Florida, Georgia, Hawaii, Idaho, Illinois, Indiana, Iowa, Kansas, Kentucky, Louisiana, Maryland, Michigan, Massachusetts, Minnesota, Mississippi, Missouri, Montana, Nebraska, Nevada, New Hampshire, North Carolina, North Dakota, Oklahoma, Oregon, Pennsylvania, South Carolina, South Dakota, Tennessee, Texas, Vermont, Washington, West Virginia, Wisconsin, Wyoming

Red Bag Only - California, Maine, New Jersey, New Mexico, New York, Ohio, Rhode Island, Virginia

Site Specific - Utah

Always wear hand and eye protection when handling liquid medical waste and follow all regulations.

Refer to the Safety Data Sheet (SDS) for detailed safety information. SDS can be printed from www.deroyal.com

For more information about Premium Solidifier+™, please call 1-800-251-9864 or visit www.deroyal.com

The DeRoyal logo, Improving Care. Improving Business, and the Premium Solidifier+ logo are trademarks or registered trademarks of DeRoyal Industries, Inc.

DeRoyal Industries, Inc.
 200 DeBusk Lane | Powell, TN 37849
 1.800.DEROYAL or 865.938.7828
www.deroyal.com | Rev. 9/18 Reprint #0-2366



Improving Care. Improving Business.®